

Refer a Patient

Phone: 833-570-3829
Fax: 512-676-5350



Patient Transfer System

Our Patient Transfer System (PTS) offers a multitude of capabilities within a single product. Designed to allow a single caregiver to easily and safely transport and reposition a patient, the chairs can perform the functions of stretcher, mobile chair, and custom seating.

H-250

- Weight capacity of 250 lbs
- Provides safe transfer to and from bed
- Single caregiver transfers made simple



Seat width: 17", Overall width: 24"

I-Series

- Weight capacity of 400 or 700 lbs
- Extended width to comfortably accommodate larger patients
- Limitless number of positions



I-400 chair: Seat width: 21.5", Overall width: 29.5"
I-700 chair: Seat width: 27", Overall width: 35"

Watch our chair videos:



Transfers are as easy as 1-2-3-4!
<http://www.humancare.se/us/product-category/products-convertiblechairs/>



See what others have to say about our Convertible Chairs.
<http://www.humancare.se/us/product-videos-transportolar/>

Who Qualifies For The Patient Transfer System?

Qualifications

In order to qualify for a Patient Transfer System, a patient must meet the following requirements.

- Patient requires a transfer between bed and either a chair, wheelchair or commode
- Without the use of a sling lift, patient is bed-confined
- Patient requires supine/lateral transfers
- Patient must be non-weight bearing

Common Qualifying Diagnoses

- Alzheimers/Dementia
- Amputations (AKA, BKA)
- Multiple sclerosis (MS)
- Joint Contractures
- Decubitus Ulcers (Stage II or greater on trunk or pelvic region)
- Degenerative Joint Disease (DJD)
- Osteoarthritis
- Paralysis (Paraplegia, Quadriplegia)
- Rheumatoid Arthritis
- Seizure Disorder
- Stroke

* Please note these are just examples of the common diagnoses. There are more that will qualify.



Medicare Covered and Coded
E1035 and E1036




1 Place the chair by bed in stretcher position.




2 Attach bed sheet to transfer bar.



3 Transfer patient by turning handle.



4 Incline chair and you're ready to go!



Benefits of the Patient Transfer System

The work required of caregivers has become increasingly dangerous due to the constant demand to move and transfer dependent patients. With that in mind, we are determined to:

- Reduce preventable incidents
- Provide the safest transfer possible
- Help you to enhance the quality of care you can provide
- Create a more positive patient experience
- Make your job easier

” *The convertible chair allows one single caregiver to perform safe patient transfers helping reduce the risk of injury for both patient and caregiver.* ”



99946-Tri-fold-brochure-Rev02 [2019-02]



HUMAN CARE
Your Life. Your Way.

Distributed by:

HUMAN CARE DIRECT

8006 Cameron Street, Suite K
Austin, TX 78754
Phone: 512-476-7199
Fax: 512-676-5350
Email: info.us@humancaregroup.com

Patient Transfer System